ONE MORE STEP IN ULTRASONIC SURGERY...

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SWEEP-TORSIONAL MODE 3D

Software unico al mondo

MEDICAL DEVICES AND ADVANCED TECHNIQUES

TIPS

*The only ones with over 100 T-Black tips
Skilled tips for specific procedures*

SURFACES TREATMENTS

T-BLACK hand touch Butterfly Effect

T-COR new surface abrasive finishing process. We are eliminating diamond from tips

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The continuous commitment of Esacrom Research and Development Department led to today realization of more than 100 specific tips for each surgery procedure, all T-Black finished.

Every tip is an example of the added value and effort that every day Esacrom makes to provide innovative solutions to surgical procedures problems.

T-BLACK

ESACROM
ELECTROMEDICAL AND DENTICAL SOLUTIONS

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*The one and only with over 100
T-BLACK tips*

Micro-affilatura al posto della diamantatura
*Micro-sharpening instead of diamond
coating*

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
*BUTTERFLY EFFECT: wetting the tip with
two drops of water, we give you the emotion
of touching on a butterfly wings!*

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“Comparison analysis on T-Black inserts have found out a reduction on the overheating of the tissue and an increased resistance at wearing out and corrosion”.

Prof. A. Scarano, Department of Oral Science, Nano e Biotechnologies University “G. D’Annunzio” Chieti (Italy).

CARATTERISTICHE
FEATURES

ELEVATA DUREZZA
 HIGHER HARDNESS

RESISTENZA ALL’ATTACCO CHIMICO
 BETTER RESISTENCE TO CHEMICAL
 ATTACK

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SUPERFICIE LISCIA
 SLIDE SURFACE

ANTIRIFLESSO
 ANTIREFLECTIVE SURFACE

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T-COR is a special treatment which performs for materials subtraction and not the addiction. Sending surface advantages are really remarkable (no residual are released).

The special surface micro-abrasive technique is further improved thanks to T-BLACK treatment, which makes it even more slippery.

ES020XLT
L'UNICA CHE
PERMETTE 3 AZIONI
1 UNICA PUNTA:
fresare, limare ed abrader.
ES020XLT
THE ONLY ONE ALLOWING
3 ACTIONS WITH 1 TIP
milling, filing and abrading.

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
The advantages:

- High efficiency;
- Operating accuracy;
- Time reduction.

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*Brand new “USER FRIENDLY”
SURGERY TECHNIQUES used for the
differentiated implant site preparation
SUS-SURGERY ULTRASONIC SITE*

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Sp: THICKNESS: it indicates the thickness of cutting tips;

Ø: DIAMETER: it indicates the diameter in the spherical and conical shaped tips;

LO: OPERATIVE LENGHT: lenght of the working area;

LA: LENGHT FROM THE CORNER: it indicates the measure of the tip from the edge to the top of the tip;

TP: DEPTH MARKINGS: laser marking of depth markings in mm from the top of the tip;*

PM: LASER MARKING: cue ball laser marking in mm from the end of the tip;*

U: SUGGESTED POWER;

V: SUGGESTED VIBRA;

P: SUGGESTED WATER PUMP;

MAX POWER: maximum power which can be used with the insert. At maximum power tip's breakeage risks is higher.

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

Our tips are sold together with our devices, in kit or in individually. The label attached on the packaging shows quantity, code and production batch.

*LE PUNTE POSSONO AVERE LE MARCATURE LASER IDENTIFICATE CON TACCHE O PALLINI. TIPS MAY HAVE LASER MARKINGS IDENTIFIED BY NOTCHES OR CUE BALLS.
E' OBBLIGATORIO CONTROLLARE LE MISURAZIONI IT'S MANDATORY THE MEASUREMENTS CHECKING

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
Trapezium shaped for remodelling and biopsy. SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	70

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Flat scalpel frontly sharpened for osteotomies, bone remodelling. SHARPENED and SAW SHAPED.

U (1)	50
V (2)	80
P (3)	100
MAX POWER (4)	70

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Flat scalpel for osteotomies in thin crest and cut teeth apices Blocs harvesting. SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	70

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
5 teeth saw for high efficiency osteotomies for biopsy. SAW SHAPED.

U (1)	45
V (2)	80
P (3)	100
MAX POWER (4)	70

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
Flat scalpel, sharpened on three edges, used for osteotomies and extraction of included teeth. SHARPENED.

U (1)	45
V (2)	80
P (3)	100
MAX POWER (4)	70

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Round scalpel for remodelling and biopsy. SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	70

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Flat scalpel for osteotomy in delicate structures and sinus lift surgery. DIAMOND COATED 150 micron.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Spherical shaped tip for delimitation of antrostomy. DIAMOND COATED 150 micron.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	60

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
*Spherical shaped tip for delimitation of antrostomy.
DIAMOND COATED 150 micron.*

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Conical shaped tip for the mobilization of the bone in antrostomy UNSHARPENED.

U (1)	20
V (2)	40
P (3)	100
MAX POWER (4)	50

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
Round scalpel for the initial unstick of sinus membrane UNSHARPENED.

U (1)	12
V (2)	20
P (3)	100
MAX POWER (4)	20

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Round scalpel for the final unstick of sinus membrane UNSHARPENED.

U (1)	12
V (2)	20
P (3)	100
MAX POWER (4)	20

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Conical shaped tip for the mobilization of the bone in antrostomy, inner neck diamonded. DIAMOND COATED 150 micron.

U (1)	20
V (2)	40
P (3)	100
MAX POWER (4)	50

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ES004BT


smooths the edges and mobilizes the bone block 2 ACTIONS, 1 TIP

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Reducing the risk of membrane's perforation thanks to less sharp edges.

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ES005T
great sinus lift with side access.

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Flat scalpel, right bended sharpened on three edges, used for osteotomies and extraction of included teeth. SHARPENED.

U (1)	30
V (2)	80
P (3)	100
MAX POWER (4)	40

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Flat scalpel, sharpened on three edges, left bended, used for osteotomies and extraction of included teeth. SHARPENED.

U (1)	30
V (2)	80
P (3)	100
MAX POWER (4)	40

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
Flat scalpel, sharpened on three edges, used for osteotomies and extraction included teeth. SHARPENED.

U (1)	45
V (2)	80
P (3)	100
MAX POWER (4)	70

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
Pointed shaped tip for simple and complex extractions. UNSHARPENED.

U (1)	20
V (2)	80
P (3)	100
MAX POWER (4)	25

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
Pointed shapened tip for simple and complex extractions. MICRO-SHARPENED.

U (1)	20
V (2)	80
P (3)	100
MAX POWER (4)	25

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Triangular flat shaped tip, sharpened on two edges for high efficiency osteotomies, extraction, teeth cutting. SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Conic shaped tip, micro-sharpened with 12 teeth, used for extraction. MICRO-SHARPENED.


U (1)	35
V (2)	80
P (3)	100
MAX POWER (4)	60

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Pointed shaped tip used for the initial preparation of the implant site. DIAMOND COATED 30 micron.

U (1)	20
V (2)	80
P (3)	100
MAX POWER (4)	25

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
Pointed shaped tip used for the initial preparation of the site. MICRO-SHARPENED.

U (1)	20
V (2)	80
P (3)	100
MAX POWER (4)	25

INITIAL PREPARATION OF THE IMPLANT SITE

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Pointed shaped tip used for the drilling of the seat before positioning of blades.

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Conic shaped tip, micro-sharpened with 8 teeth, used for the initial preparation of the implant site.

U (1)	35
V (2)	80
P (3)	100
MAX POWER (4)	60

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Pointed shaped tip Used for initial implant site preparation. UNSHARPENED with special surface T-COR treatment.

U (1)	20
V (2)	80
P (3)	100
MAX POWER (4)	25

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SUS: SURGERY ULTRASONIC SITE

Ultrasonic site differentiated preparation technique biologically guided conceived and achieved by Esacrom Research & Development Department.

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

Conical tip for implant site preparation. MICRO-SHARPENED.

U (1)	40
V (2)	90
P (3)	100
MAX POWER (4)	60

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
Conical shaped tip for implant site preparation. MICRO-SHARPENED.

U (1)	40
V (2)	90
P (3)	100
MAX POWER (4)	60

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
Conical shaped tip for implant site preparation and FLARING MICRO-SHARPENED.

U (1)	40
V (2)	90
P (3)	100
MAX POWER (4)	60

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Conical shaped tip for implant site preparation, MICRO-SHARPENED.

U (1)	40
V (2)	90
P (3)	100
MAX POWER (4)	60

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
Conical shaped tip for implant site preparation and FLARING. MICRO-SHARPENED.

U (1)	40
V (2)	90
P (3)	100
MAX POWER (4)	60

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
INNOVATION

The Research and Development Department Esacrom, after performing the technological upgrading and softwares on the SURGYSONIC console family, introducing the exclusive "SWEEP-TORSIONAL MODE 3D" software not mechanical, unique in the world, designs, implements and creates the 1st complete T-Black kit for the implant site preparation "SUS-SURGERY ULTRASONIC SITE", in cooperation of:

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
Conic shaped tip, micro-sharpened with 8 teeth, used for the initial preparation of the implant site. MICRO-SHARPENED.

U (1)	35
V (2)	80
P (3)	100
MAX POWER (4)	60

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Crown used for the implant site preparation. MICRO-SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Conical shaped tip, micro-sharpened with 8 teeth, used for the initial preparation of the implant site. MICRO-SHARPENED.

U (1)	35
V (2)	80
P (3)	100
MAX POWER (4)	60

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Conical shaped tip for the preparation of the implant site. MICRO-SHARPENED.

U (1)	35
V (2)	80
P (3)	100
MAX POWER (4)	40

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*Conical shaped tip for the preparation of the implant site.
MICRO-SHARPENED.*

U (1)	35
V (2)	80
P (3)	100
MAX POWER (4)	40

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*Conic shaped tip for the preparation of the implant site.
MICRO-SHARPENED.*

U (1)	35
V (2)	80
P (3)	100
MAX POWER (4)	40

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NO DIAMOND YES MICRO-SHARPENED

High efficiency, operating accuracy, reduction of working, less bone overheating.

FINAL PREPARATION OF THE IMPLANT SITE

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Olive shaped tip used for the preparation of the last two millimeters of the implant site, near to the sinus or to the nerve. DIAMOND COATED 40 micron.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Cylindrical shaped tip, used for the preparation of the last two millimeters of the implant site, near to the sinus or to the nerve. DIAMOND COATED 150 micron.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Spherical shaped tip used for the preparation of the last two millimeters of the implant site, close to the maxillary sinus or to the nerve DIAMOND COATED 150 micron.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	60

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Cylindrical shaped tip, used for the preparation of the last two millimeters of the implant site near to the sinus or to the nerve. MICRO-SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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FOR CORTICOTOMIES

Orthodontic purpose by mini-invasive access to be used also by left-handed operators.

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Sharpened flat tip on a right sided bend SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Sharpened flat tip on a left sided bend SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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
Open bilaterally bite, in first class with occlusal dental crowding and contraction of the transvers diameter of the superior maxilla.

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

Treated with orthodontic therapy combined with corticotomies to obtain an expansion of about 4 mm. in the premolar region and adequate alveolar morphology .

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Half-moon shaped tip for corticotomy of the interdicular hypoplastic septa. SHARPENED.

U (*)	40
V (*)	80
P (3)	100
MAX POWER (4)	70

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Half-moon shaped left bended tip for corticotomy of the interdicular hypoplastic septa. SHARPENED.

U 1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Half-moon shaped right bended tip for corticotomy of the interradiolar hypoplastic septa.

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Flat scalpel sharpened on three edges for corticotomy subapical o in the upper zygomatic-maxillary area with cortical > 3mm. SHARPENED.

U (1)	45
V (2)	80
P (3)	100
MAX POWER (4)	70

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Flat scalpel sharpened on three edges for corticotomy subapicalo in the upper zygomatic-maxillary area with cortical > 3mm. SHARPENED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Saw shaped 4 teeth, used for superficial decortication.

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Crown shaped tip for core drilling SAW SHAPED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Left bended tip used for the preparation of the seat for the retrograde filling in frontal elements.

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Right bended tip used for the preparation of the seat for the retrograde filling in frontal elements.

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**FOR ROOT CANAL PREPARATION AND
BROKEN TOOLS REMOVAL**

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*Pointed shaped tip, used for
root canal preparation and
broken tools removal.
UNSHARPENED.*

U (1)	4
V (2)	0
P (3)	100
MAX POWER (4)	6

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Pointed shaped tip, used for root canal preparation and broken tools removal. DIAMOND COATED 40 micron.

U (1)	4
V (2)	0
P (3)	100
MAX POWER (4)	6

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
Pointed shaped tip, used for root canal preparation and broken tools removal. UNSHARPENED.

U (1)	4
V (2)	0
P (3)	100
MAX POWER (4)	6

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
Used for root canal preparation / endofile holder 0,8 mm. UNSHARPENED.

U (1)	4
V (2)	0
P (3)	100
MAX POWER (4)	6

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Used for the cleaning of the fork. Left bended. UNSHARPENED.

U (1)	8
V (2)	0
P (3)	100
MAX POWER (4)	10

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Used for the cleaning of the fork. Right bended. UNSHARPENED.

U (1)	8
V (2)	0
P (3)	100
MAX POWER (4)	10

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Pointed shaped tip for radicular sanding in non surgical periodontal therapy. With special T-COR surface treatment.

U (1)	4
V (2)	0
P (3)	100
MAX POWER (4)	15

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Pointed shaped elongated tip for radicular sanding in non surgical periodontal therapy. With special T-COR surface treatment.

U (1)	4
V (2)	0
P (3)	100
MAX POWER (4)	10

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T-COR is a special treatment which performs materials subtraction and not the addition. Sanding surface advantages are really remarkable (no residual are released).

The special surface micro-abrasive technique is further improved thanks to T-BLACK treatment, which makes the external even more slippery.

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LESS FRICTION MORE SANDING

thanks to the special T-Black finishing of the tips

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
*Used for upper gengiva scaling.
UNSHARPENED.*

U (1)	8
V (2)	0
P (3)	100
MAX POWER (4)	10

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Used for lower gingiva scaling. UNSHARPENED.

U (1)	8
V (2)	0
P (3)	100
MAX POWER (4)	10

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
Used for upper gingiva scaling. UNSHARPENED.

U (1)	8
V (2)	0
P (3)	100
MAX POWER (4)	10

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Used for gingiva scaling. UNSHARPENED.

U (1)	8
V (2)	0
P (3)	100
MAX POWER (4)	12

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
*Conical shaped tip used for the finishing of the abutment.
DIAMOND COATED 30 micron.*

U (1)	15
V (2)	0
P (3)	100
MAX POWER (4)	20

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*Conical shaped tip used for the finishing of the abutment.
DIAMOND COATED 70 micron.*

U (1)	15
V (2)	0
P (3)	100
MAX POWER (4)	20

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*Half lateral sphere used in the
interproximal preparation.
DIAMOND COATED 30 micron.*

U (1)	20
V (2)	0
P (3)	100
MAX POWER (4)	25

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*Half lower sphere used in the
interproximal preparation.
DIAMOND COATED 30 micron.*

U (1)	20
V (2)	0
P (3)	100
MAX POWER (4)	25

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*Flat scalpel diamond coated for
interproximal preparation. DIA-
MOND COATED 30 micron.*

U (1)	20
V (2)	0
P (3)	100
MAX POWER (4)	25

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
*Pointed shaped Tip used for
interproximal preparation.
DIAMOND COATED 30 micron.*

U (1)	20
V (2)	0
P (3)	100
MAX POWER (4)	25

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Saw shaped tip with 5 teeth, used for osteotomies at a high efficiency. SAW SHAPED.

U (1)	45
V (2)	80
P (3)	100
MAX POWER (4)	70

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Saw shaped tip with 4 teeth, used for osteotomies at a high efficiency. SAW SHAPED.

U (1)	45
V (2)	80
P (3)	100
MAX POWER (4)	70

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Saw shaped with 4 teeth left bended, used for osteotomies at a high efficiency. SAW SHAPED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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*Saw shaped tip with 4 teeth right bended, used for osteotomies at a high efficiency. **SAW SHAPED.***

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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*Flat scalpel tip sharpened on three edges, used for osteotomies at high efficiency. **SHARPENED.***

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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*Flat scalpel left bended tip sharpened on three edges, used for osteotomies at high efficiency. **SHARPENED.***

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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*Flat scalpel right bended, sharpened on three edges, used for osteotomies at high efficiency. **SHARPENED.***

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Flat scalpel, saw shaped on three edges, used for high efficiency osteotomies.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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Flat scalpel, sharpened on three edges, used for high efficiency osteotomies.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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
5 teeth saw for high efficiency osteotomies for biopsy SAW SHAPED.

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	70

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
*Triangular shaped tip, Sharpened on two edges for high efficiency osteotomies, extraction, teeth cutting. **SHARPENED.***

U (1)	40
V (2)	80
P (3)	100
MAX POWER (4)	50

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The One & Only Ultrasonic inserts kit concepted for Zygomatic implants

ES052XZT	ES052XZLT	ES3.5ZT	ES3.5ZLT
L. 74 mm	L. 94 mm	L. 74 mm	L. 94 mm
ø 3,0 mm	ø 3,0 mm	ø 3,5 mm	ø 3,5 mm


Power (U)	50
Vibration (V)	80
Max Power* (U)	60

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Conical shaped tip, 8 teeth micro-sharpened. For first hole. MICRO-SHARPENED

U (1)	50
V (2)	80
P (3)	100
MAX POWER (4)	60

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truncated conical shaped tip, 8 teeth microsharpened. For final hole. MICRO-SHARPENED

U (1)	50
V (2)	80
P (3)	100
MAX POWER (4)	60

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



L'UNICA LUCE A FASCIO DIFFUSO **THE ONE AND ONLY WITH** **SPREAD LED LIGHT BEAM**

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Best worldwide Led handpiece, manufactured with incomparable led light intensity, easily manageable and without shadings.

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User friendly upgrading software, to use the maximum advantages of brand new techniques and surgery procedures Simply use a USB flash drive and connect it to the port on the back of the SURGYSONIC medical device

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surgysonic MOTO
BREAKING NEW GROUNDS

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Medical Device for ultrasonic surgery ready for micromotor upgrade.

surgysonic MOTOS

GROW WITH TECHNOLOGY
INVESTING IN TIME

5 Funzioni programmabili: Surgery, Normal, Sweep, Turbo, Endo.
5 functions: Surgery, Normal, Sweep, Turbo, Endo.

SURGYSONIC MOTO

ESACROM
electronics and medical devices

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SURGYSONIC MOTO is the new approach to hard tissues surgery. Thanks to the combination into a single device of the "classic" high speed micromotor and the modern "ultrasonic" technologies, Surgysonic Moto can be used not only in the DENTAL field, but also in GENERAL SURGERY (eg. maxillo-facial surgery, neurosurgery, etc.).

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During the hole creation, at the achieving of the pre-setted torque, the motor makes 1 autoreverse turn and re-start again in the drilling direction (Function adjustable via software).

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surgysonic II
BREAKING NEW GROUNDS

NON SOLO CHIRURGIA OSSEA **BONE SURGERY AND MUCH MORE**

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With surgysonic II a new generation of surgical ultrasonic devices was born. A brand new version, more powerful in performance and with a new user- friendly software.

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surgysonic II duo
TWO IS BETTER THAN ONE

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*Basic version with one handpiece, upgradable
with the second one.*

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With BMS IMPLANT II a new generation of surgical brushless micromotor was born. A brand new version more powerful in performance and with a user friendly software. Just one contra-angle for implantology and for endodontics applications with reduction 20:1 velocity from 10 up to 2000 rpm constant torque from 1 ncm up to 77 ncm. Visualization at the end of the work of the average torque, the seconds of work and number of blocks of the instrument in endo function.

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
"RECIPROCATING" TECHNIQUE

Cw/Ccw: function forwards-backwards programmable by software a step of 12°.

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The tips showed in this brochure can be bought either in convenient kit and individually.

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Thanks to new Esacrom peristaltic pump the insertion of the irrigation tube is easier

SURGYMAN L: manipoled ad ultrasuoni
SURGYMAN L: ultrasonic led handpiece

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NEW DESIGN:

- compact, ergonomic, damage - resistant;
- new finishing techniques: anti-scratches, unfadeable colours;
- possible choice between 6 different colours.

COLORS AND DESIGN

Esacrom and Sardi are partners in the development, innovation and design of our complete product lines. Sardi is the leader in the "innovation business" supporting Esacrom's staff in the continuous development of new products and technologies for the international medical-surgical field.

DESIGN By Enrique Luis Sardi & Esacrom Lab Team

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ERGONOMIC DESIGN: to avoid the stress of the foot.

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
NO TOUCH: working without touching the housing of the device.

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RHEOSTAT FUNCTION: possible customized calibration of the footswitch, based on the foot pressure.

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HALL EFFECT: functioning based on magnetic field; no electronic parts undergo wearing out.

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SURGYMAN: ULTRASONIC HANDPIECE FOR BONE SURGERY (3)

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SURGYMAN L: SULTRASONIC LED HAND PIECE FORE BONE SURGERY. (3)

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IM300MAII: MICROMOTOR HANDPIECE FOR IMPLANTOLOGY no contrangle

SURGYSUPP: SUPPORTO MANIPOLO DESTRO E SINISTRO

SURGYSUPP: HANDPIECE HOLDER, RIGHT AND LEFT

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AST SAC: HOLDER FOR PHISIOLOGIC LIQUID

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*Power source: 230Vac 50/60Hz – 115V.50 Hz
Power consumption Nominal: 170 VA
Length of the handpiece cable 2000 mm
Capacity of the hydraulic circuit: 0 ml/min 70ml/min
Pneumatic footswitch
Dimensions 300x280x130 mm
Console weight: 4.5 kg
Ultrasonic data:
Vibration 20/200 micron
Functions: Surgery, Normal, Endo, Turbo, Sweep
Program: 10 each function
Frequency: 22000-35000 hz
Power 50w (70 w with turbo function and sweep function)
Weight of the ultrasonic handpiece including cable and connector: 182 gr.
Micromotor data:
Brushless motor with hall sensors
Speed: 200-40.000 rpm
Torque: 77 Ncm
Ncm from 0,3 to 77
Functions: drill, ream, implantology, Prothesis, Endo
Programms. 10 each function
Weigth of the micromotor handpiece: 350 gr
Possible exclusion of the pump
Reverse
Autoreverse
Waterproof case for transport
Dimensions of the packed parcel (cardboard): cm 24,5x44x61,5*

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*Power source: 230Vac 50/60Hz – 115V 50 Hz
Power consumption Nominal: 170 VA
Length of the handpiece cable 2000 mm
Capacity of the hydraulic circuit: 0 ml/min 70ml/min
Pneumatic Footswitch
Dimensions 300x280x130 mm
Console weight: 4,5 kg
Ultrasonic data:
Vibration 20/200 micron
Functions: Surgery, Normal, Endo, Turbo, Sweep
Program: 10 each function
Frequency: 22000-35000 hz
Power 50w (70 w with turbo function and sweep function)
Weight of the ultrasonic handpiece including cable and connector: 182 gr.
Micromotor data (in the SURGYSONIC MOTOS version is optional):
Brushless motor with hall sensors
Speed: 200-40.000 rpm
Torque: 77 Ncm
Ncm from 0,3 to 77
Functions: Drill, Ream, Impl., Prosthesis, Endo
Programms: 10 each function
Weigth of the micromotor handpiece: 350 gr
Possible exclusion of the pump
Reverse
Autoreverse
Waterproof case for transport
Dimensions of the packed parcel (cardboard: cm 24,5x44x61,5).*

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Power consumption Nominal: 170 VA

Length of the handpiece cable 2000 mm

Capacity of the hydraulic circuit: 0 ml/min 70ml/min

Pneumatic footswitch

Dimensions 300x280x130 mm

Console weight: 4,5 kg

Micromotor data:

Brushless motor with hall sensors

Speed: 200-40.000 rpm

Torque: 77Ncm

Ncm from 0,3 to 77

Functions: drill, ream, implantology, Prothesis, Endo

Programms. 10 each function

Weigth of the micromotor handpiece: 350 gr

Possible exclusion of the pump

Reverse

Autoreverse

Tissue case for transport

Dimensions of the packed parcel (cardboard): cm 24,5x44x61,5

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Muletto in sostituzione entro 48 ore.
Substitutive device for temporary replacement during repair.

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Sconto su accessori e consumabili.
Special discount on the purchase of accessorise and consumables.

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ESACROM BELIEVES IN THE IMPORTANCE OF CONTINUING TRAINING AND FOR THIS PURPOSE ORGANIZES DIFFERENT TYPES OF COURSES IN ITALY AND AROUND THE WORLD.

For Proper inserts' use, the parameters managing, knowledge of operating mode, acquisition of Innovative techniques and procedure for the implant site preparation is suggested the joining to Esacrom dedicated training courses.

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BASIC *to whom who approaches ultrasonic surgery. For small groups, in the company or in Clinic, to learn the usage of ultrasonic devices, and its physical and mechanical application.*

ADVANCED *for individual or for small groups, possibility of customization. Study on the several surgical application through videos, Live Surgery, hands-on on simulators, surgery on patients with the support of a tutor.*

PHYSIC LAB *theoretical-practical study on the main physical or mechanical principal of ultrasound and on the techniques of usage of surgical tips based on the different materials, finishing and shapes.*

ANATOMIC DISSECTION *anatomic theoretical-practical study. Hands-on with our devices are made on cadavers and in comparison with traditional techniques.*

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At the end of the course all participants will receive a Certificate of Attendance

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“Esacrom has been founded in 1999 thanks to the sharing of ideas, about the scientific research and new technologies development, by some technical experts in electronics.

The solutions proposed so far have proved often of excellence in the various sectors in which we operate.

The continuous research of our engineers has produced over time numerous and tangible benefits for professionals and patients.

The improved operational procedures combined the high efficiency of appliances and tools have allowed simplification of surgery procedures and on the other a significant improvement in life patients’ quality.

Esacrom a step forward in ultrasonic surgery.”

p. 16/17 ENTRY LEVEL

..... ES001T
 ES005T
 ES006T
 ES007T
 ES009T
 ES010T

p. 18/21 SINUS LIFT

..... ES002T
 ES008AT
 ES015T
 ES004T
 ES003AT
 ES003BT
 ES004BT
 ES005T

p. 22/25 EXTRACTION

..... ES009T
 ES009RT
 ES009LT
 ES009XT
 ES012T
 ES012XT
 ES009NT
 ES052XFT

**p. 26/27 INITIAL PREPARATION
 OF THE IMPLANT SITE**

..... ES012AT
 ES012XT
 ES012GNT
 ES052XGT
 ES012CT
 ES030ACT

**p. 28/37 SUS-SURGERY
 ULTRASONIC SITE**

..... ES009NT
 ES052XGT
 ES02.8T
 ES03.2T
 ES03.6T
 ES0SV1T
 ES04.0T
 ES04.4T
 ES0SV2T

p. 38/41 CILINDRIC DRILLING

..... ES052XGT
 ES040T
 ES041T
 ES043T
 ES044T
 ES045T
 ES046T

..... ES047T
 ES048T
 ES049T

p. 42/43 CONICAL DRILLING

..... ES052XGT
 ES056T
 ES057T
 ES058T
 ES059T
 ES060T

**p. 44/45 FINAL PREPARATION
 OF THE IMPLANT SITE**

..... ES08AT
 ES08BT
 ES020T
 ES020XT
 ES020XLT

p. 46/47 ESA INCISA

..... ES007DRT
 ES007DLT

p. 48/49 ESA CORTIX

..... ES007GT
 ES007GLT
 ES07GRT
 ES009T
 ES091T
 ES007LT

p. 50/51 **BONE SCRAPING & PICKING**

- ES0034T
- ES0035T
- ES07ST
- ES010T
- ES001T

p. 52/53 **ENDO APICECTOMY**

- ES033LT
- ES033RT
- ES031T
- ES032T
- ES002T
- ES015AT

p. 54/55 **ENDODONTICS**

- ES023T
- ES023DT
- ES024T
- ES024DT
- ES011T

p. 56/58 **PERIODONTOLOGY**

- ES003DT
- ES030LT
- ES030RT
- ES030LDT
- ES030RDT

- ES021T
- ES012T
- ES012CT
- ES030ACT

P. 59/61 **SCALER**

- ES016T
- ES018T
- ES019T
- ES021T
- ES012T

P. 62/63 **ABUTMENT FINISHING**

- ES050FT
- ES050GT
- ES051FT
- ES051GT
- ES052FT
- ES052GT

P. 64/65 **CAVITY**

- ES08C1T
- ES08C2T
- ES08C3T
- ES08C4T
- ES09DT
- ES012AT

P. 66/69 **MAXILLO**

- ES071T
- ES007AT
- ES007LT
- ES007RT
- ES090T
- ES091T
- ES092T
- ES009ST
- ES009L5T
- ES007L5T
- ES009SFT
- ES009NT

P. 70/73 **ZYGOMA ULTRASONIC DRILLING**

- ES052XZT
- ES052XZLT
- ES3.5ZT
- ES3.5ZLT

3	R&D innovations	64	Cavity
4	Tips	66	Maxillo
5	T-Black	70	Zygoma ultrasonic drilling
8	Micro-sharpening	75	Guide lines
9	T-Cor	76	Hi-Led
10	New techniques	77	Software upgrade
11	Led Light Handpiece	78	Surgysonic MOTO
12	Sweet-torsional mode 3D	79	Surgysonic MOTOS
14	Info	80	Surgysonic MOTO
16	Entry level	82	Surgysonic II
18	Sinus lift	85	Surgysonic II DUO
22	Extraction	86	Surgysonic II DUOS
26	Initial preparation of the implant site	87	BMS Implant II
28	SUS-Surgery Ultrasonic Site	90	Kit & medical devices
38	Cilindric drilling	92	Design & colors
42	Conical drilling	93	Footswitch
44	Final preparation of the implant site	94	Equipment & accessories
46	Esalncisa	96	Series Equipment
48	EsaCortix	98	Technical data
50	Bone scraping & picking	102	Customer care
52	Endo apicectomy	103	Training
54	Endodontics	104	Special thanks
56	Periodontology	105	Mission
59	Scaler	108	Tips index
62	Abutment finishing	110	Index

Ringraziamo Thanks to:

Design: ESACROM LAB TEAM and Enrique Luis Sardi

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Text and communication: ESACROM LAB TEAM

ESACROM[®]
electronics and medical devices

Via Zambrini 6/A - 40026 Imola (Bo) ITALY
Tel. +39 0542 643527 - Fax +39 0542 482007
esacrom@esacrom.com
www.esacrom.com

Surgysonic[®] name is registered in china:

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